



## SWITCH CONTROL

### How to Access:

1. Open SETTINGS
2. Tap GENERAL
3. Tap ACCESSIBILITY
4. Under 'INTERACTION', tap 'SWITCH CONTROL'.
5. Turning on SWITCH CONTROL will allow you to customize various features, such as: Number of switches connected, scanning style and type, **the interaction of the switch(es), and any sound effects or speech.**
  - a. You can add a new switch via Bluetooth, connecting it with a port, or through head movements using the front facing camera.
  - b. Once a switch is added or selected, there are multiple ways to use: manual selection, item scanning, or point scanning.
  - c. **Once a specific tab/ feature is selected, a menu selection will open depending on what feature was selected.**

### How to Activate:

1. Ensure the 'SWITCH' tab is 'ON'
2. Review all features and ensure that customizable items are set up to the child's preferences.
3. Begin using SWITCH ACCESS!

### Videos Tutorials:

<https://www.youtube.com/watch?v=TH540kuBDwo>

<https://www.youtube.com/watch?v=GQKEE9n1lk>

